

Product Datasheet

HEA125 antibody (orb27221)

Description

Mouse Monoclonal antibody againsts HEA125


Species/Host


Mouse

Conjugation

Unconjugated

Tested

IHC-Fr, IHC-P

Applications

Immunogen

Living human colon carcinoma cell line HT-29

Immunohistochemical staining of prostate...

Preservatives

Antibody solution in stabilizing phosphate buffer pH 7.3. Contains 0.09 % sodium azide** and rAlbumin. The volume is sufficient for at least 50 immunohistochemical tests (100 µl working solution / test). Use appropriate antibody diluent e.g. BIOLOGO Art. No. PU002, if further dilution is required.

Storage

Storage: 2-8°C

Note

For research use only

Application notes

Positive Control: Appendix. Pre-Treatment: Pre-treatment with pronase 0.1% or unmasking fluid C (Art. No. DE000). Incubation Time: 60 min at RT

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID

P16422

Hazard
Information

*These antibodies are intended for in vitro research use only. They must not be used for clinical diagnostics and not for in vivo experiments in humans or animals. ** The preservative sodium azide is known to be poisonous and potentially hazardous to health. It should be handled only by trained staff. Despite of the product's low azide concentration it must be handled with care. Dispose according to regional rules!

Dilution Range

(RTU) neat

Expiration Date

12 months from date of receipt.