

Product Datasheet

**HPV (Type 1, 6, 11, 16, 18, 31) antibody
(orb27223)**

Description

Mouse Monoclonal to HPV (Type 1, 6, 11,

Species/Host

Mouse

Conjugation

Unconjugated

Tested

IHC-Fr, IHC-P

Applications
Immunogen

SDS-disrupted bovine papillomavirus type 1

Preservatives

Antibody solution in stabilizing phosphate buffer pH 7.3. Contains 0.09 % sodium azide**. The volume is sufficient for at least 30 immunohistochemical tests (100 µl working solution / test). Use appropriate antibody diluent e.g. BIOLOGO Art. No. PU002 if required.

Storage

Storage: 2-8°C

Note

For research use only

Application notes

Positive Control: Infected cervical biopsy. Pre-Treatment: temperature antigen retrieval using 0.01 M citrate retrieval solution (pH 6.0) is recommended. Incubation Time: 60 min at RT

Isotype

IgG

Clonality

Monoclonal

Hazard
Information

*These antibodies are intended for in vitro research use only. They must not be used for clinical diagnostics and not for in vivo experiments in humans or animals. ** The preservative sodium azide is known to be poisonous and potentially hazardous to health. It should be handled only by trained staff. Despite of the product's low azide concentration it must be handled with care. Dispose according to regional rules!

Dilution Range

(RTU) neat

Expiration Date

12 months from date of receipt.