

Product Datasheet

HEA125 antibody (orb27259)

Description

Mouse Monoclonal antibody againsts HEA125


Species/Host

Mouse

Conjugation

Unconjugated

Tested

FC, IHC-Fr, IHC-P

Applications
Immunogen

Living human colon carcinoma cell line HT-29

Preservatives

Antibody solution in stabilizing phosphate buffer pH 7.3. Contains 0.09 % sodium azide**. The volume is sufficient for at least 100 immunohistochemical tests (100 µl working solution / test). Use appropriate antibody diluent e.g. BIOLOGO Art .No. PU002.

Concentration

50 µg/ml

Storage

Storage: 2-8°C

Note

For research use only

Application notes

Positive Control: Appendix. Pre-Treatment: Pre-treatment with pronase 0.1% or unmasking fluid C (Art. No. DE000). Incubation Time: 60 min at RT

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID

P16422

Hazard
Information

*These antibodies are intended for in vitro research use only. They must not be used for clinical diagnostics and not for in vivo experiments in humans or animals. ** The preservative sodium azide is known to be poisonous and potentially hazardous to health. It should be handled only by trained staff. Despite of the product's low azide concentration it must be handled with care. Dispose according to regional rules!

Dilution Range

(liquid conc.) 1:10-1:50

Expiration Date

12 months from date of receipt.

Immunohistochemical staining of Prostate...