

Product Datasheet

IL18bp antibody (HRP) (orb27453)

Description	IL18bp antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Interleukin-18-binding protein (29-193AA)
Target	IL18bp
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Z0M9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.