

Product Datasheet

CML (N-Epsilon)-Carboxymethyl-Lysine antibody (orb27499)

Description

Mouse Monoclonal to CML (N-Epsilon)-

Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	N-Epsilon-(Carboxymethyl)-Lysine (CML)-HSA
Preservatives	Protein G affinity purified antibody from ascites in stabilized buffer, containing 50% Block Ace? (Casein-containing solution, Dainippon Co.) and 0.1% ProClin? (Rohm & Haas) as a preservative
Concentration	0.25 mg/ml
Storage	Storage: -20°C
Note	For research use only
Application notes	Incubation Time: 60 min at RT or 18 hr at 2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Hazard Information	These antibodies are intended for in vitro research use only. They must not be used for clinical diagnostics and not for in vivo experiments in humans or animals.
Dilution Range	(liquid conc.) ELISA 0.1-0.5 µg/ml; IHC 2 µg/ml
Expiration Date	12 months from date of receipt.