

Product Datasheet

Vitamin D3 Receptor antibody [Out of stock] (orb27513)

Description	Rat monoclonal antibody to Vitamin D3 Receptor
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	FA, IHC-P, IP, WB
Immunogen	Partially purified chicken vitamin D receptor
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Dilution Range	(liquid conc.) 1:500-1:1000 (IHC), 1:1000-1:2000 (WB)
Expiration Date	12 months from date of receipt.