

Product Datasheet

Platelet factor 4 PF4 Antibody (orb27644)

Description

Platelet factor 4 PF4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human PF4 (85-101aa DLQAPLYKKI IKLLES).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

10845 MW

Uniprot ID
P02776
Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By HeatWestern blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



— 40KD



— 20KD




— 14KD



— 8KD

Western blot analysis of Recombinant Hum...



5ng(Lane2) and 2.5ng(Lane3) using CXCL4 ...