

Product Datasheet

Goat Anti-Rabbit IgG (5nm) antibody (Gold) [Out of Stock] [Out of stock] (orb27691)

Description	Goat Anti-Rabbit IgG (5nm) antibody
Species/Host	Goat
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	Rabbit IgG (whole molecule)
Form/Appearance	Concentrated, Liquid
Storage	At 4°C for one year.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Electron Microscopy 1:20-50 Immunohistochemistry(Paraffin-embedded Section) 1:50-100 Immunohistochemistry(Frozen Section) 1:50-100 Immunocytochemistry 1:50-100
Expiration Date	12 months from date of receipt.