

Product Datasheet

Goat Anti-Rabbit IgG (10nm) antibody (Gold) [Out of Stock] [Out of stock] (orb27692)

Description	Goat Anti-rabbit IgG (H+L) secondary antibody, 10 nm Colloidal Gold Conjugate, for detection, locali...
Species/Host	Goat
Reactivity	Rabbit
Tested Applications	EM, ICC, IF, IHC-Fr, IHC-P
Immunogen	Rabbit IgG (whole molecule)
Form/Appearance	Concentrated, Liquid
Storage	At 4°C for one year.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Electron Microscopy 1:20-50 Immunohistochemistry(Paraffin-embedded Section) 1:50-100 Immunohistochemistry(Frozen Section) 1:50-100 Immunocytochemistry 1:50-100
Expiration Date	12 months from date of receipt.