



Product Datasheet

Cy3 Conjugated Avidin Antibody (orb27703)

Description	Cy3 Conjugated Avidin Antibody
Conjugation	Cy3
Tested Applications	FC, IF
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



IF analysis of human intestinal cancer t...



IF:Beta Catenin was detected in paraffin...