

Product Datasheet

Goat Anti-Mouse IgG (H+L) Secondary Antibody, Cy3 Conjugated (orb27704)

Description

Goat Anti-Mouse IgG (H+L) Secondary



Species/Host

Goat

Conjugation

Cy3

Tested

FC, IF

Applications

Concentration

1mg/ml

Storage

At -20°C for one year from date of receipt.
Avoid repeated freezing and thawing. Protect from light.

Note


For research use only

Clonality


Polyclonal

Expiration Date


12 months from date of receipt.



IF:PCNA
was
detected in
human
intestinal...



IF:PCK was
detected in
human
intestinal ...



Beta
Tubulin was
detected in
an
immunocy...