



Product Datasheet

Goat Anti-Rabbit IgG (H+L) Secondary Antibody, Cy3 Conjugated (orb27705)

Description	Goat Anti-Rabbit IgG (H+L) Secondary Antibody, Cy3 Conjugated
Species/Host	Goat
Conjugation	Cy3
Tested Applications	FC, IF
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



IF:Collagen IV was detected in rat kidney...



IF:MCM2 was detected in human intestinal...