

Product Datasheet

THOC1 antibody (orb27796)

Description	Mouse monoclonal antibody to THOC1 is expressed in various cancer cell lines. It participates in an ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	THOC1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG3,k
Clonality	Monoclonal
MW	75666
Uniprot ID	Q96FV9
NCBI	NP_005122.2
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

