

Product Datasheet

S100A8 antibody (Biotin) (orb27807)

Description	S100A8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	H. pylori
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Helicobacter pylori Cytotoxicity-associated immunodominant antigen protein (918-1147AA)
Target	cagA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P80200
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.