

Product Datasheet

SOX9 antibody (orb27816)

Description	Mouse monoclonal antibody.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SOX9
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM,k
Clonality	Monoclonal
MW	56137
Uniprot ID	P48436
NCBI	NP_000337.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

