


Product Datasheet

Ribosome-inactivating protein cucurmosin antibody (orb27831)

Description	Ribosome-inactivating protein cucurmosin antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Cucurbita moschata Ribosome-inactivating protein cucurmosin protein (1-244AA)
Target	Ribosome-inactivating protein cucurmosin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84531
Dilution Range	WB:1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Cucurbita mosch...