

Product Datasheet

C1QL2 antibody (orb2784)

Description	C1QL2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human C1QL2/C1QTNF10 (101-187/187 aa)
Target	C1QTNF10
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27 kDa
Uniprot ID	Q7Z5L3
Expiration Date	12 months from date of receipt.