

Product Datasheet

Ribosome-inactivating protein charybdin antibody (HRP) (orb27849)

Description	Ribosome-inactivating protein charybdin antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Drimia maritima Ribosome-inactivating protein charybdin protein (1-257AA)
Target	Ribosome-inactivating protein charybdin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84786
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.