




Product Datasheet

PRMT2 antibody (orb27874)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PRMT2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49042
Uniprot ID	P55345
NCBI	NP_001229794.1 , NP_001526.2 , NP_996845.1 , NP_001229793.1
Dilution Range	IHC-P: 1:50-100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): 293 c...



Immunohistochemical analysis of formalin...