

Product Datasheet

Pka-C1 antibody (orb27908)

Description	Pka-C1 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Drosophila melanogaster cAMP-dependent protein kinase catalytic subunit 1 protein (2-353AA)
Target	Pka-C1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12370
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.