

Product Datasheet

Protein X antibody (orb27932)

Description	Protein X antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Woodchuck hepatitis B virus Protein X protein (1-141AA)
Target	X
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17401
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.