




Product Datasheet

BCA3 antibody (orb27946)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AKIP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23114
Uniprot ID	Q9NQ31
NCBI	NP_001193577.1 , NP_001193575.1 , NP_001193576.1 , NP_065693.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of BCA3 antibody i...



Western blot analysis of BCA3 antibody i...