

Product Datasheet

Ribosome-inactivating protein karasurin-C antibody (orb27973)

Description	Ribosome-inactivating protein karasurin-C antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Trichosanthes kirilowii Ribosome-inactivating protein karasurin-C protein (22-270AA)
Target	Ribosome-inactivating protein karasurin-C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24478
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.