

Product Datasheet

AAP1 antibody (orb27989)

Description	AAP1 antibody
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Saccharomyces cerevisiae Alanine/arginine aminopeptidase protein (1-389AA)
Target	AAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37898
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.