

Product Datasheet

Nuclease P1 antibody (orb28007)

Description	Nuclease P1 antibody
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Penicillium citrinum Nuclease P1 protein (1-270AA)
Target	Nuclease P1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24289
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.