




Product Datasheet

ATP5H antibody (orb28047)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ATP5PD
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18491
Uniprot ID	O75947
NCBI	NP_006347.1 , NP_001003785.1
Dilution Range	IHC-P: 1:50-100, FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
ATP5H
antibo...



Western blot
analysis of
ATP5H
antibody ...