

Product Datasheet

phyA antibody (Biotin) (orb28083)

Description	phyA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Aspergillus niger 3-phytase A protein (24-467AA)
Target	phyA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P34752
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.