

Product Datasheet

PPDK antibody (orb28101)

Description	PPDK antibody
Species/Host	Rabbit
Reactivity	Parasite
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Entamoeba histolytica Pyruvate, phosphate dikinase protein (1-342AA)
Target	PPDK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37213
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.