

Product Datasheet

Dnmt3A/Dnmt3A2 antibody (orb28111)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DNMT3A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	101858
Uniprot ID	Q9Y6K1
NCBI	NP_783329.1 , NP_715640.2 , NP_072046.2 , NP_783328.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

