

Product Datasheet

JMJD2D antibody (orb28167)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KDM4D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	58603
Uniprot ID	Q6B016
NCBI	NP_060509.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

