

Product Datasheet

PMAP36 antibody (Biotin) (orb28174)

Description	PMAP36 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Pig Antibacterial peptide PMAP-36 protein (130-166AA)
Target	PMAP36
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49931
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.