

Product Datasheet

vacA antibody (FITC) (orb28182)

Description	vacA antibody (FITC)
Species/Host	Rabbit
Reactivity	H. pylori
Conjugation	FITC
Immunogen	Recombinant Helicobacter pylori Vacuolating cytotoxin autotransporter protein (37-245AA)
Target	vacA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55981
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.