

Product Datasheet

POLR2J antibody (Biotin) (orb28188)

Description	POLR2J antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human DNA-directed RNA polymerase II subunit RPB11-a protein (1-117AA)
Target	POLR2J
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P52435
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.