

Product Datasheet

Delta-ctenitoxin-Pn1a antibody (HRP) (orb28192)

Description	Delta-ctenitoxin-Pn1a antibody (HRP)
Species/Host	Rabbit
Reactivity	Other
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Phoneutria nigriventer Delta-ctenitoxin-Pn1a protein (35-82AA)
Target	Delta-ctenitoxin-Pn1a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P59368
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.