




Product Datasheet

SLC5A12 antibody (orb28215)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	SLC5A12
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67647
Uniprot ID	Q1EHB4
NCBI	NP_848593.2
Dilution Range	WB: 1:2000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Flow cytometric analysis of CEM cells showing increased fluorescence in the ME population.