

Product Datasheet

fkpA antibody (Biotin) (orb28222)

Description	fkpA antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli O6:H1 FKBP-type peptidyl-prolyl cis-trans isomerase FkpA protein (26-270AA)
Target	fkpA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P65764
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.