

Product Datasheet

UL48 antibody (FITC) (orb28245)

Description	UL48 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Human herpesvirus 1 Tegument protein VP16 protein (1-490AA)
Target	UL48
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04486
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.