

Product Datasheet

Histidine-rich protein PFHRP-II antibody (HRP) (orb28260)

Description	Histidine-rich protein PFHRP-II antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Plasmodium falciparum Histidine-rich protein PFHRP-II protein (31-289AA)
Target	Histidine-rich protein PFHRP-II
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05227
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.