

Product Datasheet

SES1 antibody (Biotin) (orb28276)

Description	SES1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Saccharomyces cerevisiae Serine--tRNA ligase, cytoplasmic protein (1-462AA)
Target	SES1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07284
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.