

Product Datasheet

mutT antibody (HRP) (orb28288)

Description	mutT antibody (HRP)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli 8-oxo-dGTP diphosphatase protein (1-129AA)
Target	mutT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08337
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.