

Product Datasheet

mpt63 antibody (Biotin) (orb28323)

Description	mpt63 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA)
Target	mpt63
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WIP1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.