

Product Datasheet

PHTF2 antibody (orb28337)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target PHTF2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 88757

Uniprot ID [Q8N3S3](#)

NCBI [NP_001120832.1](#), [NP_001120830.1](#), [NP_065165.3](#), [NP_001120829.1](#), [NP_001120831.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

