

Product Datasheet

lepB antibody (orb28383)

Description	lepB antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Signal peptidase I protein (78-324AA)
Target	lepB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00803
Dilution Range	ELISA : 1: 5000 - 10000
Expiration Date	12 months from date of receipt.