

Product Datasheet

Histone H3-MeK9 antibody (orb28415)

Description Rabbit polyclonal antibody to histone H3-MeK9 belongs to the histone H3 family that constitutes the ...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target H3-3A

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 15328

Uniprot ID [P84243](#)

NCBI [NP_005315.1](#), [NP_002098.1](#)

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

