


Product Datasheet

ox LDL antibody (orb2844)

Description	ox LDL antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	Mouse Cu ²⁺ oxidised low density lipoprotein
Target	ox-LDL
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31 kDa
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis
of HepG2 cell
using ox-LDL

...