

Product Datasheet

gshA antibody (orb28467)

Description	gshA antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Glutamate--cysteine ligase protein (1-518AA)
Target	gshA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A6W9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.