


Product Datasheet

ZFYVE28 antibody (orb28487)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	ZFYVE28
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	96490
Uniprot ID	Q9HCC9
NCBI	NP_001166130.1 , NP_001166127.1 , NP_066023.2 , NP_001166129.1 , NP_001166131.1 , NP_001166128.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...