



Product Datasheet

S100Z antibody (orb28537)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	S100Z
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	11620
Uniprot ID	Q8WXG8
NCBI	NP_570128.2
Dilution Range	FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB453, K562...