




Product Datasheet

NGRN antibody (orb28547)

Description	Rabbit polyclonal antibody to WDR32, which is also named as DCAF10 interacts with DDB1. WDR32 could ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	NGRN
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	32408
Uniprot ID	Q9NPE2
NCBI	NP_001028260.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cell line



Flow cytometric analysis of A549 cells (...)