

Product Datasheet

OR10X1 antibody (orb28549)

Description	Rabbit polyclonal antibody to CD109 is a GPI-anchored member of the alpha-2-macroglobulin (A2M) and ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OR10X1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36446
Uniprot ID	Q8NGY0
NCBI	NP_001004477.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

